

Do you know about the Oasis barcode system?

By Tracey Penn

The Oasis barcoding system is an add-on component to the Oasis Enterprise system. The components include a Sony Barcode Label Printer, a Barcode Reader and the Oasis Barcoding software.

Autoclave Records

Oasis streamlines the recording of autoclave data.


Set-up and printing



Oasis prints-out barcode labels which are then affixed to instrument packs as they come out of the autoclave. Up to 99 different label formats can be printed - ideal for different autoclave machines. Two labels are printed for each barcode number and the number is unique for that machine. One label goes in your manual book and the other onto the pack. Facilities are also provided on each barcode label to define your own headings - i.e. manual recording of date/time, cycle number and operator etc.

Usage

As you open packs in the surgery, use the barcode reader to read the barcode label. Oasis will add the autoclave record as a special clinical note with the usual auditing protection. This means you have a permanent computer record for that autoclave pack against each patient's records viewable on screen within the Clinical Notes and Charting windows (With our optional Query Report Writer you can actually generate reports on autoclave records by date range).



Stock control records

Oasis Barcoding also works with stock control records. Using the normal Oasis invoicing screen, stock items can be read in using the barcode reader.

Oasis recognizes the normal supplier bar codes on your sales items and enables you to link the stock item to an Oasis item code. This means faster and more accurate data entry procedures.

Stock items without supplier barcodes can also be processed by Oasis giving you the ability to manually key in the bar code and Oasis printing the bar code labels for you.

Patient labels

Because the Sony Label Printer uses regular single sticky labels in addition to printing barcode labels, it can also be used for many other label printing functions.

The Barcode Kit also includes a range of preset label formats for use in printing patient address labels, referring doctor labels, Teledex labels and the like. We also give you the tools to design your own labels!

For more information on Oasis, please call 1300-889-668 or in New Zealand call 0800-772-216.